



Transport System Folded

The Porter Transport System is a fast, simple securing system with 2 never at wheels and eight position handles, offering an all purpose transport solution. Greater leverage and easier carrying reduces fatigue and chance of injury.

Features

- Universal design, compatible with any type and style of rigid litter
- Eight position adjustable handles
- Fast ratchet securing system locks basket stretcher with 2 straps and 2 rubberized steel hooks/strap to cinch it down
- Quick release pins on each side of the handle locks them into position with ease; Each pin has a wire lanyard
- Design allows the wheels to be quickly and easily removed for compact storage and transport
- Two securing straps (2 hooks per strap) are each rated to 4.5 kN (1,000 lbs.) and use a strong and durable Ratchet-Locking device easily operated in all weather conditions

Specifications

Porter Transport System		
Length	37.5 in	95 cm
Width	28 in	71 cm
Depth	19.5 in	50 cm
Weight	41 lbs	18.5 kg
Rated Load Limit	500 lbs	227 kg