


Gazelle Series Rescue Litters

11-0105, Regular
11-0106, Tapered
11-0107, Confined Space



The **Gazelle Rescue Litter Series** is intended for basic rescue and first response. Available in 3 models: Regular, Tapered, and Confined Space.

Features

- Tubular 1.9 cm top frame and lower rails
- HDPE - High density polyethylene netting
- Black powder coat finish
- Reinforced head and foot ends
- Rectangular Model -- Allows any backboard
- Class 1 medical device

Specifications

Gazelle Series Rescue Litter		
Length (ALL)	82.5 in	209.6 cm
Width (REGULAR, TAPERED)	24.5 in	62.2 cm
Width (CONFINED)	18.3 in	46.4 cm
Depth (ALL)	7 in	17.8 cm
Weight (REGULAR, TAPERED)	33 lb	15 kg
Weight (CONFINED)	31 lb	14 kg
Rated Load Limit (ALL)	900 lb	4 kN