

Pinnacle Rescue Litter

11-4001, Split Apart
11-4000, Regular



The **Pinnacle Rescue Litter** is the world's **lightest rigid stretcher** for confined space rescue. Its narrow design allows for entry into and egress from most confined spaces. Available as One-Piece or Split-Apart.

Features

- StratLoad™ feature for vertical lift
- HDPE - High Density Polyethylene netting
- UHMW-PE Back Support - No more splinters from a wooden back support, and it is extremely durable
- All TIG welded frame
- Titanium material is extremely light weight
- Confined Design
- Class 1 medical device

Specifications

Pinnacle Rescue Litter		
Length	82.5 in	209.6 cm
Width	19 in	48.3 cm
Width - Split Apart	19.5 in	49.5 cm
Depth	7 in	17.8 cm
Depth - Split Apart	7.3 in	18.4 cm
Weight	16 lbs	7 kg
Weight - Split Apart	17 lbs	8 kg
Rated Load Limit	900 lbs	408 kg