



Hard Eye Litter Bridle  
15-0105



Soft Eye Litter Bridle  
15-0106

For use in rescue situations where raising and lowering of a basket stretcher is required during technical rescue operations. Use of Mil-Spec hardware throughout reduces wear on the web.

## Features

- Fully adjustable litter bridles for use with any style litter
- Steel "D" ring on the end of each leg of the hard eye litter bridle
- Webbing loop on the end of each leg of the soft eye litter bridle
- Parachute buckles used on each of the four legs allow them to be adjusted independently -- stretcher can easily be set at the best angle for patient comfort
- Includes mesh carrying bag
- Carabiners sold separately
- Class 1 medical device

## Specifications

Traverse Adjustable Litter Bridles		
Weight (HARD EYE)	3.4 lb	1.5 kg
Weight (SOFT EYE)	2.8 lb	1.3 kg
Rated Load Limit (ALL)	5000 lb	22 kN