

SAM
MEDICAL

BLEEDING CONTROL KIT

BE PREPARED. PRESERVE LIFE.

Severe blood loss, often caused by motor vehicle accidents, acts of violence and natural disasters, can result in death within minutes. The SAM Medical Bleeding Control Kit includes essential medical devices designed to limit blood loss in emergency and traumatic situations. After calling emergency personnel, this life-preserving kit may help stabilize the injured in such events.



FEATURED TOURNIQUET

SAM XT EXTREMITY TOURNIQUET is recommended by the U.S. Department of Defense's Committee on Tactical Combat Casualty Care (CoTCCC).

Engineered for rapid application, SAM XT's innovative design normalizes the number of windlass turns needed compared to traditional tourniquets, enabling an easier, faster training of the product and intuitive use.



KIT CONTENTS:

1 SAM XT EXTREMITY TOURNIQUET

Slack is the main cause of failed tourniquet application. SAM XT's TRUFORCE™ Buckle technology auto-locks at a predetermined amount of circumferential force, eliminating nearly all tourniquet slack. Simply click, twist, and secure for a life-preserving tourniquet application. The SAM XT is on the list of approved articles by the DoD for the Stop The Bleed Initiative.

2 MEDICAL GLOVES

Nitrile protective medical gloves - Latex free, powder free, size L.

3 MEDICAL SHEARS

Stainless steel medical shears allow the first responder to expose the wound and provide for appropriate treatment more efficiently.

4 COMPRESSED GAUZE

Designed to treat wounds involving moderate to severe blood loss. Can be secured in place with Emergency Bandage. 100% cotton, 6-ply, sterile.

5 EMERGENCY BANDAGE

Pressure bandage can be used directly over a wound or to hold gauze dressing in place.

6 INSTRUCTION CARD


Includes detailed instructions for use and a diagram for proper tourniquet placement.

7 PERMANENT MARKER


To indicate time of tourniquet application on patient.



CONTENTS



SAM XT
Extremity Tourniquet



Medical
Gloves

Medical
Shears



Compressed
Gauze



Permanent
Marker



Emergency Bandage
6 in wide (15.2 cm)



Instruction Card

PRODUCT SPECIFICATIONS

PACKAGED SIZE & WEIGHT*

Length: 10.1 in (25.8 cm)
Width: 7 in (18 cm)
Depth: 3.2 in (8 cm)
Weight: 10.9 oz (310 g)

OPERATING SPECIFICATIONS

After calling emergency personnel, this life-preserving kit may help stabilize the injured in traumatic events.

AVAILABLE MODELS

Vacuum packed (KT901-V-EN)
Non-Vacuum packed (KT901-NV-EN)
Bleeding Control Pack -contains x4
Bleeding Control Kits (KT920-M4-EN)

GLOBAL SALES

U.S. Toll Free: 800.818.4726
U.S. Local: 503.639.5474
U.S. Fax: 503.639.5425
sammedical.com
sammedical.store



SAM Medical Products®
12200 SW Tualatin Rd., Ste 200
Tualatin, OR 97062 USA



* Size and weight vary by configuration and manufacturing process.
SAM® is a registered trademark in the United States and other countries.
©2020 SAM Medical Products. All rights reserved.