

# SAE Packages

**Make Your Ambulance Compliant** with the SAE Ambulance Safety Requirements that went into effect July 1, 2015 with **Ferno's SAE Compliance Packages.**



**iNjX™** Integrated Patient Transport & Loading System™ with iNjLine Fastening System with **integrated charging** (includes IV Pole & Sof-Net).



**POWERFlex+®** Powered Cot with **Stat Trac®** Cot Fastening System and **ICS Fastener Charging System** (includes IV Pole & Sof-Net).



**35X PROFlexx®** Manual X Frame Cot with **Stat Trac®** Cot Fastening System (includes IV Pole & Sof-Net).

Each of our SAE Compliance Packages includes a cot, fastener, and charging system.\*

After more than four years of testing and research, the Society of Automotive Engineers (SAE) has published their "Recommended Practice" concerning the forces that **an ambulance cot, fastening system and patient restraints should be able to withstand in a 30 mile-per-hour crash.** These requirements will be adopted into new standards, including GSA KKK-A-1822F and NFPA-1917. **Traditional "Antler and Rail" fastening systems will NOT meet these new requirements.**

Ferno compliance packages, featuring the iNjLine Fastening System and the Stat Trac® Cot Fastening system, have been dynamically tested and are **fully compliant** with the new requirements from SAE. **Upgrade your ambulance now to protect your crews and patients from injury.**

**CALL 877.733.0911 FOR PRICING**

# SAE Packages

The Ferno® iNjLine Fastening System for iNjX, and the Stat Trac® Cot Fastening System have been dynamically tested and are **fully compliant with new ambulance safety requirements from the Society of Automotive Engineers (SAE).**



## iNjLine Fastening System for iNjX

- In-line design that eliminates drift for easier loading and unloading
- iNjX can be loaded from any angle or on an incline, without any lifting
- Integrated charging for iNjX.



## Stat Trac Cot Fastening System

- Straight in-track design eliminates cot drift and the need to steer the cot in, minimizing strain on angled surfaces
- A single-hand release and wide-angle entrance guide make cot loading quick & easy
- Can accept the Ferno POWERFlexx+ST Powered Ambulance Cot
- Can accept 35XST PROFlexx Manual X-Frame Cot
- Optional ICS Integrated In-the-Fastener Charging



\*ICS is an optional accessory for POWERFlexx

**FOR MORE INFORMATION CONTACT:**