



Features

- Designed for use with rescue litters
- 65 lb. of buoyancy to keep patient afloat
- Conveniently stores inside rescue litter
- Compatible with Titan Series, Gazelle Series, Advantage and Model 71 Basket Stretchers
- PVC (vinyl) coated polyester - latex free



Flotation Collar, 2 Piece

The lightweight Traverse **Flotation Collar, 2 Piece** provides flotation assistance during patient handling in water.

It's flexible, 2 piece design allows for quick and easy installation on a wide variety of rescue litters.

Specifications

Flotation Collar, 2 Piece		
Length	51 in	130 cm
Width	4 in	10 cm
Depth	8 in	20 cm
Weight	3 lbs	1.4 kg
Buoyancy	65 lbs	29.4 kg