



The **Titan Rescue Litter** is a complete rethinking of the classic Stokes military-style litter.

The top rail is 1" (2.5 cm) diameter, making the stretcher easier to handle when carried for extended periods. Class 1 medical device.

Features

- Available in Stainless Steel
- Wider design for bariatric patients
- Compatible with any backboard
- Patented StratLoad™ attachment points for secure attachment of lifting bridles
- Replaceable HDPE netting allows water and air flow
- Designed to surpass the US Military specification MIL-L-37957A. Tested to and beyond a 1,134 kg (2500 lb) static load test threshold
- Large 1" top tube for easy handling



StratLoad™ feature for vertical lift

Specifications

Titan 32 SS (Stainless Steel)		
Length	83.5 in	212.1 cm
Width	32.4 in	82.2 cm
Depth	7.4 in	18.7 cm
Weight	42 lbs	19 kg
Rated Load Limit	2500 lbs	11 kN