



Features

- Titanium frame, lighter for bariatric use
- Patented StratLoad™ attachment points for secure attachment of lifting bridles
- Replaceable HDPE netting allows water and air flow plus a durable UHMW-PE back support
- Split-apart feature with LocSafe™ connectors provide fail-safe connection and allows the stretcher to be assembled in about 15 seconds
- Designed to surpass the US Military specification MIL-L-37957A. Tested to and beyond a 1,134 kg (2500 lb) static load test threshold
- Large 1" top tube for easy handling
- Tapered design enables easier access in confined areas



StratLoad™ feature for vertical lift

The Tirol Titan TI Rescue Litter is a Titanium, TIG-welded frame construction Rescue Litter, featuring our patented LocSafe Connection system and StratLoad™ attachment points.

Class 1 medical device.

Specifications

Tirol Titan TI, With Plates		
Length	83.5 in	212.1 cm
Width	23.3 in	59.1 cm
Depth	7.6 in	19.4 cm
Weight	21 lbs	9.5 kg
Rated Load Limit	2500 lbs	11 kN

Weighed with Accessories