



The **Bradco Rescue Litter System** is ideally suited for **industrial applications** of any kind. Its built-in restraint system makes it easy to use and it is so robust that it will take the rigors of any environment. The patient can quickly be restrained for either a vertical or a horizontal evacuation.

Features

- Well-suited for any type of rescue, including industrial applications
- Designed for horizontal and vertical use
- Heavy-duty, all-aluminum frame protects patient and has a 2,500-lb strength rating
- Double-flap cover allows patient to be transported securely in almost any position
- Polyethylene shell is resistant to chemicals and glides smoothly over various types of terrain
- Suspended restraints system allows for fast immobilization and comfortable evacuation

Specifications

Bradco Rescue Litter System		
Length	80.4 in	204.2 cm
Width	24.4 in	61.9 cm
Depth	8.6 in	21.9 cm
Weight	44 lbs	20 kg
Rated Load Limit	2500 lbs	11 kN