



Mule II Litter Wheel



Securing Straps



Compact

The Mule II Litter Wheel with an all terrain wheel, fast, simple securing system and 8 position handles offers an all purpose transport solution. Greater leverage and easier carrying reduces fatigue and chance of injury.

Features

- Universal design, compatible with any type and style of rigid litter
- 8 position adjustable handles
- Single-wheel design ideal for narrow trails, allows for agility, versatility and ease of movement in the roughest terrain including bogs, ruts, roots
- Large, lightweight ATV wheel offers greater tracking ability and shock absorption for smooth and comfortable patient transport
- Fast ratchet securing system locks basket stretcher with 2 straps and 2 rubberized steel hooks/strap to cinch it
- New quick release pins on each side of the handle locks them into place

Specifications

11-2210 11-2210G

Mule II (Deployed)

Length	33 in	84 cm
Width	26 in	66 cm
Depth	29 in	74 cm
Weight	36 lbs	16 kg

11-2220 11-2220G

Mule II with Handles (Deployed)

Length	96 in	244 cm
Length (Stored)	26 in	66 cm
Width	27.5 in	70 cm
Width (Stored)	32 in	81 cm
Depth	29 in	74 cm
Depth (Stored), (No Handles)	19 in	48 in
Weight	46lbs	21kg
Rated Load Limit (No Handles)	500 lbs	227 kg