

# Fluid Resistant Restraints

*Ferno Essentials*

## Spec Sheet



417 Series



Model 773



Model 430P



Model 430

## Features

- 100% Polyvinyl Chloride (PVC) webbing
- Polyester scrim 17 x 18 Thread
- Tear Strength 140 lbs
- Hemmed Tab
- Fill Strength 120 lbs
- Fluid-Resistant Rated
- -40 degree Cold Crack
- Not crash-tested. Do not use with ambulance cots

## Specifications

|             |                                                                      |
|-------------|----------------------------------------------------------------------|
| PART NUMBER | PRO 90001 (Model 430), PRO 90004 (Model 430P), PRO 90006 (Model 773) |
|-------------|----------------------------------------------------------------------|

|        |          |      |
|--------|----------|------|
| LENGTH | 177.8 mm | 7 in |
|--------|----------|------|

|       |         |      |
|-------|---------|------|
| WIDTH | 50.8 mm | 2 in |
|-------|---------|------|

|             |                                                                      |
|-------------|----------------------------------------------------------------------|
| PART NUMBER | PRO 90002 (Model 430), PRO 90003 (Model 430P), PRO 90005 (Model 773) |
|-------------|----------------------------------------------------------------------|

|        |        |      |
|--------|--------|------|
| LENGTH | 127 mm | 5 in |
|--------|--------|------|

|       |         |      |
|-------|---------|------|
| WIDTH | 50.8 mm | 2 in |
|-------|---------|------|

|             |                        |
|-------------|------------------------|
| PART NUMBER | PRO 90007 (417 Series) |
|-------------|------------------------|

|      |          |
|------|----------|
| SIZE | Standard |
|------|----------|



Add: 2460 Tedlo Street,  
Mississauga, ON, L5A 3V3

Web: [www.ferno.ca](http://www.ferno.ca)

Ph: 800.543.3766

Fax: 905.615.1199